


Amidia presents the revolutionary

Nano Diamond Knives

Using nano diamond knives
is the most economical
way to perform a clear
corneal or a side-port
incision.

Get a diamond scalpel
at the price of a
sapphire knife.



swiss made

*sharpness
at it's peak*

amidia

Amidia LTD
Mattenstrasse 11
CH-2555 Brügg

Call +41 32 333 24 74
www.amidia.ch


Nano Diamond Knives

*sharpness
at it's peak*




drawing	product no.	description	width	length	thickness	angle
---------	-------------	-------------	-------	--------	-----------	-------

single edge




A-nano-1000	straight	1.0 mm	3.5 mm	0.15 mm	45°
A-nano-1005	straight	1.5 mm	3.5 mm	0.15 mm	45°

lance



A-nano-1100	straight	1.0 mm	3.5 mm	0.15 mm	30°
A-nano-1110	straight	1.0 mm	3.5 mm	0.15 mm	60°
A-nano-1120	straight	1.5 mm	3.5 mm	0.15 mm	70°

tri-facet




A-nano-1200	straight	1.0 mm	3.5 mm	0.15 mm	30°/30°
-------------	----------	--------	--------	---------	---------



drawing	product no.	description	width	length	thickness	angle
---------	-------------	-------------	-------	--------	-----------	-------

clear cornea



A-nano-1123A		1.8 mm	4.0 mm	0.15 mm	70°
A-nano-1127A		2.2 mm	4.0 mm	0.15 mm	70°
A-nano-1129A		2.4 mm	4.0 mm	0.15 mm	70°
A-nano-1130A		2.5 mm	4.0 mm	0.15 mm	70°
A-nano-1131A		2.6 mm	4.0 mm	0.15 mm	70°
A-nano-1133A		2.8 mm	4.0 mm	0.15 mm	70°