

REINVENT DRY EYE TREATMENT



REXONEYE

Rexon-Eye is the ultimate, **non-invasive and durable treatment** for all forms of dry eye. It works by applying low-power high-frequency electric fields, capable of stimulating the rejuvenation and natural regeneration of cells.


The principle of operation of Rexon-Eye is QMR: Quantum Molecular Resonance. This technology delivers weak alternate electric current patterns at frequencies (4 MHz and multiple) that are in resonance with tissue molecular bindings.

The Rexon-Eye therapy successfully works for both **evaporative as well as hypo-secretive dry eye**. The benefits for the patient begin just after the treatment and have been shown to be long-term. The treatment itself is **comfortable, pleasant and relaxing**.






REXONEYE



Display

Dial

Connectors



Mask

FLUIDDESIGNLAB.COM

FEATURES

- Therapy consists in 4 treatment sessions of 20 min each, at one week interval.
- Treatment is delivered by a contact electrode built in a mask, worn by patient over closed eyes.
- Instrument set-up and use are very simple and intuitive.
- Treatment is comfortable, pleasant and relaxing for the patient.
- Rexion-Eye is **CE marked since 2016** and patented in Italy (international patents pending).

BENEFITS

- The Rexion-Eye therapy successfully works for all types of dry eyes, both evaporative as well as hypo-secretive.
- Treatment results in immediate as well as long-term benefits for the patient, evidenced by subjective (OSDI questionnaire, use of tear substitutes) as well as objective measurements (e.g., TBUT, corneal staining, tear production).
- Rexion-Eye's efficacy in treating the dry eye syndrome is strongly clinically proven (Pedrotti E. et al., Br J Ophthalmol 2017).

R-Eye Rev. 01/18 – Subject to technical changes as part of product development, as well as errors and omissions.



RESONO OPHTHALMIC SRL

Headquarters
Via Leonardo da Vinci, 13
36066 Sandrigo (VI) - Italy
Mob. +39 349 - 6788770
Tel. +39 0444 - 239519

contact@resono.it
www.resono.it

Registered Office
Via Giacinto Gallina, 4
34122 Trieste - Italy
Tel. +39 040 - 764041

